

# Influx Replicate

Auf Edge RPi: docker compose exec -it influxdb /bin/bash

influx remote create \

1. -name dieters-server-ab \
2. -remote-url <https://influx.dieter.eolab.de> \
3. -remote-api-token EsB0by9yZVd0J6ShuWiloJRTjrIPCIIG-zexHNQUMHI7QlwnoZrdrEzUYG7us2HdmPemvhu01lciO4cJrl7\_eg== \
4. -remote-org-id 2d621d5bc12efe46

influx replication create \

1. -name dieters-server-DieterDataRep \
2. -remote-id 0d0874a051d59000 \
3. -local-bucket-id ecd383d225199195 \
4. -remote-bucket 22d58136410ea02f

From:

<https://wiki.eolab.de/> - HSRW EOLab Wiki

Permanent link:

<https://wiki.eolab.de/doku.php?id=eolab:dieter:influx-replicate&rev=1728985161>

Last update: **2024/10/15 11:39**

