

Automatic Speech Recognition (ASR) with vosk


Sources

vosk

- <https://alphacepei.com/vosk/>
- <https://github.com/alphacep/vosk-api>

Dataquest

- <https://github.com/dataquestio/project-walkthroughs/blob/master/microphone/microphone.ipynb>



Installation

```
conda create -n vosk python=3.9
conda activate vosk
conda install -c conda-forge jupyterlab numpy matplotlib pandas
#conda install -c conda-forge ipywidgets
#conda install -c conda-forge scipy scikit-learn
```

```
pip install vosk
pip install pyaudio
```

On Windows the vosk models are **cached here**: C:\Users\<<username>\.cache\vosk

ffmpeg

- <https://www.gyan.dev/ffmpeg/builds/>

```
import pyaudio
import wave

# Constants for audio recording
FORMAT = pyaudio.paInt16
CHANNELS = 1
RATE = 44100
CHUNK = 1024
RECORD_SECONDS = 5 # Adjust this to change the duration of the recording
OUTPUT_FILENAME = "output.wav"

def list_audio_devices():
    audio = pyaudio.PyAudio()
    devices = []

    for i in range(audio.get_device_count()):
        device_info = audio.get_device_info_by_index(i)
        devices.append(f"{i}: {device_info['name']}")

    audio.terminate()
    return devices

def get_input_device_index():
    devices = list_audio_devices()

    print("Available audio input devices:")
    for device in devices:
        print(device)

    while True:
        try:
            print("")
            print("On Becker's Dell Lat. 7330 the following works:")
            print("1: Microphone Array (Realtek(R) Au")
            print("")
            device_index = int(input("Enter the index of the desired input
device: "))
            if 0 <= device_index < len(devices):
                return device_index
            else:
                print("Invalid input. Please enter a valid device index.")
        except ValueError:
            print("Invalid input. Please enter a valid device index.")

def record_audio(device_index):
    audio = pyaudio.PyAudio()
```

```
# Open a microphone stream with the selected input device
stream = audio.open(format=FORMAT, channels=CHANNELS,
                    rate=RATE, input=True,
input_device_index=device_index,
                    frames_per_buffer=CHUNK)

print(f"Recording from:
{audio.get_device_info_by_index(device_index)['name']}")

frames = []

# Record audio in chunks and store it in frames
for _ in range(0, int(RATE / CHUNK * RECORD_SECONDS)):
    data = stream.read(CHUNK)
    frames.append(data)

print("Finished recording.")

# Stop and close the microphone stream
stream.stop_stream()
stream.close()
audio.terminate()

# Save the recorded audio to a WAV file
with wave.open(OUTPUT_FILENAME, 'wb') as wf:
    wf.setnchannels(CHANNELS)
    wf.setsampwidth(audio.get_sample_size(FORMAT))
    wf.setframerate(RATE)
    wf.writeframes(b''.join(frames))

if __name__ == "__main__":
    device_index = get_input_device_index()
    record_audio(device_index)
```

From:
<https://wiki.eolab.de/> - HSRW EOLab Wiki

Permanent link:
<https://wiki.eolab.de/doku.php?id=user:rolf001:vosk:start&rev=1694356904>

Last update: **2023/09/10 16:41**

